


Telecare Technology ...


Functionality and User Experience

Functional Components of a Telecare System



MATCH System Functional
Architecture
v5.0
Nov 2008

What Can Be sensed?



- User Behaviour
 - Activity
 - Physiological Data
 - Location
- User Metadata
 - Health Records
- User Commands
- Environmental & External Data
 - Weather Data
 - Events

User Behaviour Sensors



2 - Interconnection and Data Paths

1 - Activity Monitors

1.1 - White Goods

1.2 - Communications Systems

1.3 - Entertainment Systems

1.4 - Home Systems

1.5 - Intruder Alarm Systems

1.6 - Safety Systems

1.7 - Special Sensors

1.8 - Special Sensor Arrays

1.9 - Location Detection

2.1 - Wireless

2.2 - Wired

MATCH Floor Mat

- 2m x 1m Floor Mat to detect gait




Detecting Change

- Granularity
- Threshold
- Frequency
- Rate

Data Processing


MATCH System Functional
Architecture
v5.0
Nov 2008



Detecting Change

- What has changed?
 - What was the state before change?
 - What is the state after change?
 - Does the change matter?
- The meaning of the change is ...
 - Generic ... (high blood pressure is always worrying)
 - Person specific ... (sleeping patterns)
- The user will need a generic and a specific Profile.

User Profile Usage



Practical Exercise:

Sensing Change in your

- What changes does the client want to know about themselves?
- How will that be detected/measured
- What changes do carers need to know in order to be able to care?
- How will these things be detected/measured